INHALT

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CONTENTS

Items Supplied	9
At A Glance	9
Symbols	10
Signal Words	10
Intended Use	10
When May the Device Not Be Used?	11
Safety Notices	11
Battery	13

Attachments and Areas of Application	13
Use	13
Care and Storage	14
Troubleshooting	15
Disposal	15
Technical Data	15

ITEMS SUPPLIED

- Mini massage gun
- Fork attachment
- Flat attachment
- Cartridge attachment
- Sphere attachment
- USB cable (A-connector to USB C-connector)
- Operating instructions

Check the items supplied for completeness and the components for transport damage. If you find any damage, do not use the device but contact our customer service department.

Remove any possible films, stickers or transport protection from the device.

Never remove the rating plate and any possible warnings!

AT A GLANCE

(Picture A)

- 1 Mount for the attachment
- 2 Main unit
- 3 Handle
- 4 Intensity level control lamps
- 5  Start / Stop / Set intensity level
- 6 Charging socket
- 7 Operating lamp
- 8 Sphere attachment
- 9 Fork attachment
- 10 Flat attachment
- 11 Cartridge attachment

Dear Customer,

We are delighted that you have chosen our **MAXXMEE mini massage gun**. It operates at up to 3200 strokes per minute and thus delivers a superb powerful massage – ideal whenever you have given your all doing sport, everyday tasks or your job. Blood circulation is intensely stimulated and your flexibility is restored. You can now treat yourself to a break with a massage whenever and wherever you want because the slender and cordless power tool will fit in almost any bag. If you have any questions about the device and about spare parts/accessories, contact the customer service department via our website:

www.dspro.de/kundenservice

We hope you have a lot of fun with your **MAXXMEE mini massage gun**.

Information About the Operating Instructions

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device. The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

SYMBOLS

Danger symbols: These symbols indicate possible dangers. Read the associated safety notices carefully and follow them.

Supplementary information

Read the operating instructions before use!

Symbol for direct current

SIGNAL WORDS

Signal words indicate hazards if the associated notices are not followed.

WARNING – medium risk, may result in serious injury or death

NOTICE – may result in risk of damage to material

INTENDED USE

- The device is intended to be used for massaging the whole body and is used for relaxation. **It is not a substitute for medically required treatment.**
- The device is **not** a toy.
- The device should not be used before going to sleep. Massages have a stimulating effect, which may lead to problems getting to sleep.
- The device is for personal use only and is not intended for commercial or medical applications.
- The device should only be used as described in the operating instructions. Any other use is deemed to be improper.
- The warranty does not cover faults caused by incorrect handling, damage or attempted repairs. The same applies to normal wear and tear.

WHEN MAY THE DEVICE NOT BE USED?

- Children, people who are sleeping or unconscious and animals must not be treated with the device.
- The following people may use the device **only following consultation with a doctor:**
 - people who are fitted with a pacemaker or similar medical implants;
 - people with heart disease or who suffer from epilepsy;
 - people who have back complaints or spinal disorders;
 - people who suffer from osteomyelitis, bone inflammation (in particular if the spinal column is damaged) or from osteoporosis;
 - women who are pregnant or have just given birth;
 - people who suffer from circulation problems or thrombosis;
- people who have undergone an operation for which the treatment is not yet complete;
- people with physical abnormalities;
- people who have a malignant tumour.
- Do not use the device in the vicinity of or on:
 - injuries,
 - swelling,
 - skin irritations,
 - infected skin,
 - warts,
 - fungal infections,
 - abnormal changes in the skin,
 - bones,
 - the head,
 - the genital area,
 - body parts which have been anaesthetised or have no feeling.
- **Only trained people may use the device on joints or around joints.**

If you are unsure whether the device may be used, consult a responsible doctor.

SAFETY NOTICES

- This device is not intended to be used by people (including children) with reduced sensory or mental abilities or a lack of experience and/or knowledge, unless they are supervised by a person who is responsible for their safety or have received instruction from this person on how to use the device safely. There is a danger of injury from incorrect use!
- Do not make any modifications to the device. It may only be taken apart and/or repaired by the manufacturer, customer service department or a similarly qualified person (e.g. specialist workshop) in order to avoid any hazards.

WARNING – Danger of Injury

- Keep children and animals away from the packaging material. There is a danger of suffocation.
- Make sure that the connecting cable is always kept out of the reach of small children and animals. There is a danger of strangulation.
- If the battery leaks, avoid physical contact with the battery acid. Wear protective gloves if necessary. If you touch battery acid, rinse the affected areas with abundant clear water straight away and seek medical assistance immediately. There is a danger of burns from battery acid.
- In the interest of your own safety, inspect the device and the connecting cable for damage each time before you use it. Use the device only if it is in working order.
- Make sure that the connecting cable does not present a trip hazard.
- Do not use the device if it has malfunctioned or fallen into water. Have it checked in a specialist workshop before using it again.

WARNING – Danger from Rechargeable Batteries

- The main unit with a built-in battery must not be taken apart, thrown into a fire or exposed to high temperatures of over +40 °C. The battery could ignite or explode.

- Make sure that the device, and in particular the main unit, cannot get wet. Do not use the device if liquid has got into it. The battery in the main unit could ignite or explode.
- Do not cover the device during charging to prevent overheating.

NOTICE – Risk of Damage to Material and Property

- The device must be switched off when the device is connected to or disconnected from the mains power.
- Make sure that the connecting cable is not squashed, bent or laid over sharp edges and does not come into contact with hot surfaces.
- When the mains adapter is pulled out of the plug socket, always pull on the mains adapter and never the mains cable.
- Do not pull or carry the device by the connecting cable. Never wind the mains cable around the mains adapter.
- Use only original accessories from the manufacturer in order to guarantee that there is no interference that may prevent the device from working and to avoid possible damage.
- Protect the device from heat sources, naked flames, sub-zero temperatures, persistent moisture, wet conditions and impacts.

BATTERY

Charging the Battery

Please Note!

- A USB mains adapter corresponding to the technical data of the device is also required (see rating plate or "Technical Data"). This is not supplied with the device.
- Risk of a short circuit! Do not stick any conducting objects (e. g. a screwdriver or similar implement) into the charging socket.
- Charge the device only in closed rooms.
- Do not charge the device in rooms with high humidity.

It is essential that you charge up the battery fully before you first use the device in order to achieve the maximum charging capacity.

1. Make sure that the device is switched off. To do this, press the  button (5) for approx. 3 seconds.
2. Plug the USB-C-USB connector of the USB cable into the charging socket (6) on the device and the USB connector into the socket of a suitable USB mains adapter.
3. Connect the USB mains adapter to an easily accessible plug socket. The operating lamp (7) lights up red. When it lights up green, the battery is fully charged.
4. Pull the USB mains adapter out of the plug socket and disconnect the USB cable from the USB mains adapter and the device.

ATTACHMENTS AND AREAS OF APPLICATION

(Pictures C–F)

- **Sphere attachment (8):**
for a large muscle group such as quadriceps and gluteus maximus and for targeted relaxation for sore muscles
- **Fork attachment (9):**
for the neck, back extensors, calves or similar muscles
- **Flat attachment (10):**
suitable for relaxing the whole body
- **Cartridge attachment (11):**
for deeper tissue and for the feet, hands etc.

USE

Please Note!

- Do not use the device under a blanket or similar cover to prevent overheating.
- Do not use the device when it is connected to the mains power.
- For hygiene reasons, the sphere attachment should not be used on skin which has had cream or oil applied.
- If you experience an unpleasant feeling or tingling sensation in the hand holding the device, stop using it and consult a doctor if necessary.

1. Insert the attachment you want as far as it will go into the mount (1) on the main unit (2).
2. Press the  button (5) for approx. 3 seconds to switch on the device. The operating lamp (7) lights up green and the control lamps (4) briefly light up one after the other. The device starts to operate on intensity level 1. The operating lamp flashes green at regular intervals during use.
3. Press the  button briefly to switch between the intensity levels. The number of control lamps lit indicates which intensity level is currently set.

We recommend starting with the lowest level and gradually increasing the level of intensity.

4. Hold the device by the handle (3) and guide the attachment with gentle pressure over the area of the body you want. **Do not apply too much pressure!**

We recommend that you do not massage one area of muscle for longer than max. 15 minutes in one go.

5. To switch off the device, press the  button for approx. 3 seconds. The control lamps briefly light up one after the other and then go out. The operating lamp also goes out.

CARE AND STORAGE

Please Note!

- Do not use any sharp or abrasive cleaning agents or cleaning pads to clean the device. They may damage the surfaces.

Make sure that the device is switched off and is not connected to the mains power.

- Remove the used attachment from the main unit (2) and clean it. All attachments **apart from the sphere attachment (8)** can be cleaned with warm water and a mild cleaning agent.
- Wipe down the main unit and the sphere attachment with a slightly damp cloth.
- Store the device and all accessories in the storage case in a dry place which is protected from sunlight and is out of the reach of children and animals.

TROUBLESHOOTING

If the device does not work properly, check whether you are able to rectify the problem yourself. If the problem cannot be solved with the steps below, contact the customer service department.

Do not attempt to repair an electrical device yourself!

The operating lamp (7) does not light up red when the device is connected to the mains power.

The USB mains adapter is not plugged correctly in the plug socket or the USB connector is not plugged correctly in the charging socket (6) or in the connecting socket on the USB mains adapter.

- ▶ Correct the fit of the plug connections.

The plug socket is defective.

- ▶ Try another plug socket.

There is no mains voltage.

- ▶ Check the fuse of the mains connection.

The battery is fully charged.

- ▶ Disconnect the device from the mains power.

The device does not start.

The battery is empty.

- ▶ Charge up the battery.

DISPOSAL

Dispose of the packaging material in an environmentally friendly manner so that it can be recycled.

This device is governed by the European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). Do not dispose of the device with built-in battery as normal domestic waste, but rather in an environmentally friendly manner via an officially approved waste disposal company.

TECHNICAL DATA

Article number:	00567
ID of the operating instructions:	Z 00567 M DS V1 0820 dk
Model number:	DT-M01
Battery:	3.7 V (Li-ion battery, 3600mAh)
Power supply:	5 V, 1 A USB mains adapter (not supplied with the device)
Battery charge time:	approx. 1.5 hours
Battery operating time:	up to approx. 9 hours (with fully charged battery)

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All rights reserved.