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## LIST OF CONTENTS

Items Supplied	14
At A Glance	15
Symbols	15
Signal Words	15
Intended Use	16
When May the Device Not Be Used?	17
Safety Notices	19
Inserting the Battery	21
Use	22
Cleaning and Storage	24
Troubleshooting	24
Disposal	25
Technical Data	25

## ITEMS SUPPLIED

(Picture A)

- 1 EMS massage pen, 1 x
- 2 Tin with contact gel, 2 x
- 3 Flat attachment, 1 x
- 4 Sphere attachment, 1 x  
1 x operating instructions

Check the items supplied for completeness and the components for transport damage. If you find any damage, do not use the device but contact our customer service department.

Remove any possible films, stickers or transport protection from the device. **Never remove the rating plate and any possible warnings!**

A battery (1.5 V; AA) is also required.

---

## AT A GLANCE

(Picture B)

- 5 Tip
- 6 Pulse activity display
- 7 Intensity level display
- 8 Top button area (-): decrease intensity level, switch off
- 9 Bottom button area (+): switch on, increase intensity level
- 10 Battery compartment cover

## SYMBOLS

Danger symbols: These symbols indicate possible dangers. Read the associated safety notices carefully and follow them.

Supplementary information

Read the operating instructions before use!

Direct current

Protection class III

## SIGNAL WORDS

Signal words indicate hazards if the associated notices are not followed.

**DANGER** – high risk, results in serious injury or death

**WARNING** – medium risk, may result in serious injury or death

**NOTICE** – may result in risk of damage to material

## Dear Customer,

We are delighted that you have chosen the **MAXXMEE EMS massage pen**. It works according to the principle of EMS (Electrical Muscle Stimulation) technology: Using low current, pulses are sent to the muscles. The EMS massage pen activates rapid muscle contractions and can be used to relax the muscles. There are 9 intensity levels from which you can choose the right one that feels pleasant to you. Also ideal for carrying with you, e.g. on journeys or to training.

If you have any questions about the device and about spare parts/accessories, contact the customer service department via our website: [www.dspro.de/kundenservice](http://www.dspro.de/kundenservice)

We hope you have a lot of fun with your **MAXXMEE EMS massage pen**.

## Information About the Operating Instructions

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device. The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

## INTENDED USE

- This device is intended to be used for stimulating the skin and the tissue beneath and for electrical muscle stimulation.
- The device is only suitable for external use and not for therapeutic purposes. It is not a substitute for medical treatment. If complaints persist, consult a doctor.

- The device is for personal use only and is not intended for commercial applications.
- Use the device only for the specified purpose and as described in the operating instructions. Any other use is deemed to be improper.
- The warranty does not cover faults caused by incorrect handling, damage or attempted repairs. The same applies to normal wear and tear.

## WHEN MAY THE DEVICE NOT BE USED?

- Children, people who are sleeping or unconscious and animals must not be treated with the device.
- Do not use the device if you have a fever or illness/dental problems.
- Do not use the device while driving a vehicle.
- Do not use the device if you are attached to high-frequency surgical devices or are in the vicinity of such devices or high-frequency commercial equipment.
- Do not use the device in the immediate vicinity of devices which operate using short waves or microwaves. Maintain a minimum gap of 1 metre.
- The device must not be used near a mobile phone. Do not use a mobile phone while using the device.
- Do not use the device while you are connected to medical equipment, e.g. an ECG machine.
- Do not use the device if you have calf pain.
- Do not use the device in the first 16 weeks of pregnancy and around the belly area throughout the entire pregnancy.

Pregnant women should only use the device after consultation with their doctor.

■ The following people may use the device **only following consultation with a doctor:**

- people who are fitted with a pacemaker or similar medical implants,
- people with heart disease or who suffer from epilepsy,
- people who in the past have experienced skin irritations, itching, skin redness or other side effects from using TENS devices,
- people who suffer from circulation problems or thrombosis,
- people undergoing anticoagulant or steroid treatment or subject to special clinical conditions,
- people who have undergone an operation for which the treatment is not yet complete,
- people with a malignant tumour,
- people with nerve damage,
- people under the influence of medication,
- people who belong to a risk group or who have concerns about their health.

■ Do not use the device in the vicinity of or on:

- metallic implants,
- regions which have recently been operated on or are scarred,
- injuries and areas of swelling,
- acne, birthmarks and warts,
- injured, inflamed, (sun-)burnt, irritated or very sensitive skin,
- the rib cage close to your heart,
- the eyes, ears, thyroid glands, mucous membranes, orifices or the genital area,
- body parts with no feeling.

If you are unsure whether the device may be used, consult a responsible doctor.

## SAFETY NOTICES

- This device is not intended to be used by people (including children) with reduced sensory or mental abilities or a lack of experience and/or knowledge, unless they are supervised by a person who is responsible for their safety or have received instruction from this person on how to use the device safely. There is a danger of injury from incorrect use!
- Do not make any modifications to the device. It may only be taken apart and / or repaired by the manufacturer, customer service department or a similarly qualified person (e.g. specialist workshop) in order to avoid any hazards.

**WARNING:** Do not use this device near a bathtub, shower or any vessel containing water. Never use while bathing or showering!

### **WARNING – Danger from Batteries**

- Batteries can be deadly if swallowed. The battery and the device should therefore be kept out of the reach of children and animals. If a battery is swallowed, seek medical assistance immediately. There is a danger of burns from battery acid.
- If the battery has leaked, do not allow your skin, eyes and mucous membranes to come into contact with the battery acid. Wear protective gloves if necessary. If you touch battery acid, rinse the affected areas with abundant clear water straight away and seek medical assistance immediately. There is a danger of burns from battery acid.
- Remove a battery which has leaked immediately and clean the contacts in the device.

- Remove the battery when it is dead or if you will not be using the device for a prolonged period. This will prevent any damage caused by leaks.
- Use only the battery type specified in the technical data.
- Never immerse the device in water or other liquids and ensure that it cannot fall into water or become wet.

### WARNING – Danger of Fire

- Keep the device away from combustible materials and substances.

### WARNING – Danger of Injury

- Keep children and animals away from the packaging material. There is a danger of suffocation.
- Children must be supervised to ensure that they do not play with the device. There is a danger of injury from incorrect use.
- Do not use simultaneously with other massage or stimulation devices.
- If you experience an unpleasant feeling while using the device, stop using it immediately!
- Stop using the device if you experience skin redness or discolouration that lasts for more than 24 hours.
- Do not use the device if it has malfunctioned, fallen into water or if there is visible damage on the device or accessories. Have it checked in a specialist workshop before using it again.

## NOTICE – Risk of Damage to Material and Property

- Use and store the device only in closed rooms.
- Do not operate or store the device in rooms with high humidity.
- Do not cover the device during operation.
- Use and store the device at room temperature. Do not expose it to extreme temperatures, temperature fluctuations, fire, heat sources and direct sunlight for a prolonged period.
- Protect the device from impacts.
- Do not place heavy objects on the device or the accessories.
- In the event of damage / faults, switch the device off immediately.
- Use only original accessories from the manufacturer in order to guarantee that there is no interference that may prevent the device from working and to avoid possible damage.

## INSERTING THE BATTERY

Before you first use the pen (1), a 1.5V battery, AA (not supplied with the device) must be inserted.

Before inserting the battery, check whether the contacts in the pen and on the battery are clean, and clean them if necessary.

1. Unscrew the battery compartment cover (10).
2. Insert a 1.5V battery, AA into the pen with the positive terminal first.
3. Screw the battery compartment cover onto the pen.

Replace the battery when the power of the pen starts to wane.

## USE

---

### **Please Note!**

- The pen must always maintain skin contact with the hand that is holding it for the circuit to remain closed and the pen to work. Do not wear gloves and do not wrap anything around the pen!
  - You might also feel the electrical pulses in the hand which is holding it. However, this is harmless. If it feels unpleasant, change your grip or clasp the pen with as much of your palm as possible.
  - The attachments and the tip of the pen must not come into contact during use with metal parts, e.g. on belts, wrist watches, jewellery or similar items.
  - Do not use on wet or dirty skin or in a damp environment.
  - Do not use on skin to which cream, oil or other flammable or alcoholic products have just been applied.
  - Do not use before going to sleep. Massages have a stimulating effect, which may lead to problems getting to sleep.
  - The skin may feel warm during the treatment. This is normal and not a cause for concern.
  - We explicitly point out that with some people skin irritations or hypersensitivity may occur as a result of the electrical muscle stimulation. In this case, stop using the device immediately and consult a doctor if necessary.
- 

The EMS massage pen (1) can be used both with and without contact gel. Contact gel creates better conductivity and makes the treatment more pleasant. Without contact gel you may experience a slight burning sensation.

The pen will also work through light clothing.

**Beginners should treat one muscle group for no more than 15 minutes a day to prevent overstimulation!**

Advanced users can use the device for up to 30 minutes a day on each muscle group. We recommend no more than 2 units per day, with a rest phase of at least 4 hours between the units.

1. Screw one of the two attachments onto the tip (5) of the pen:
  - the sphere attachment (4) for targeted relaxation of fairly large tense muscle groups,
  - the flat attachment (3) for more gentle treatment of sensitive areas.

For a particularly intense massage, the pen can also be used without an attachment.

---

2. If necessary, apply a little contact gel to the area that is to be treated.
3. Take hold of the pen.
4. Press the bottom button area (9) (marked with +) to switch on the pen and increase the strength of the pulses until the display (7) shows a 2.  
The pulse activity display (6) flashes.
5. Press the tip of the pen onto the area being treated. If at first you do not feel anything, move the tip slightly until you feel a tingle. If you still do not feel a sensation, increase the pulse strength.
6. Leave the pen in position for a few minutes (for pain relief ideally 15 to 20 minutes). To treat larger areas, move the pen over these areas.
7. To reduce the pulse strength, press the top button area (8) (marked with -).
8. Switch off the pen after the treatment. To do this, keep pressing the top button area until a 0 is shown on the display. Shortly after this, the display will go out.

- Clean the attachment or the tip (see the "Cleaning and Storage" chapter).

## CLEANING AND STORAGE

---

### Please Note!

- Do not use any sharp or abrasive cleaning agents or cleaning pads to clean the device. These may damage the surfaces.
- 
- Wipe down the attachments (3, 4) and if necessary also the remaining parts of the device with a clean, damp cloth. Then dry with a soft cloth.
  - Store the device and accessories in a dry place that is protected from dust, dirt and direct sunlight and cannot be accessed by children or animals.

## TROUBLESHOOTING

If the device does not work properly, check whether you are able to rectify the problem yourself. If the problem cannot be solved with the steps below, contact the customer service department.

**Do not attempt to repair an electrical device yourself!**

The device does not work. The displays (6, 7) remain dark.

The battery is too weak or inserted incorrectly.

- ▶ Insert a new battery or insert the battery the other way round.

The stimulation is too weak or cannot be felt even though a high pulse strength has been set.

The circuit is broken.

- ▶ Hold the pen (1) with your bare hand. There must not be any hairs, clothing, jewellery etc. between the pen and the skin.

The battery is too weak.

- ▶ Insert a new battery.

The stimulation is too painful even though a low pulse strength has been set.

- ▶ Use (more) contact gel.

## DISPOSAL

Dispose of the packaging material in an environmentally friendly manner so that it can be recycled.

The battery should be removed prior to disposal of the device and disposed of separately from the device at a suitable collection point.

This device is governed by the European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). Do not dispose of this device as normal domestic waste, but rather in an environmentally friendly manner via an officially approved waste disposal company.

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All rights reserved.